

Product Datasheet

Mouse FABP3 ELISA Kit (orb439425)

Description	Mouse Fatty Acid Binding Protein 3, Muscle And Heart (FABP3) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P11404
Sensitivity	26.3pg/mL
Entrez	14077
Expiration Date	Please enquire.