

Product Datasheet

Rat FABP3 ELISA Kit (orb439427)

Description	Rat Fatty Acid Binding Protein 3, Muscle And Heart (FABP3) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P07483
Sensitivity	0.065ng/mL
Entrez	79131
Expiration Date	Please enquire.