

Product Datasheet

Human FADS3 ELISA Kit (orb439441)

Description	Human Fatty Acid Desaturase 3 (FADS3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Y5Q0
Sensitivity	0.082ng/mL
Entrez	3995
Expiration Date	Please enquire.