

Product Datasheet

Human FAIM3 ELISA Kit (orb439442)

Description	Human Fas Apoptotic Inhibitory Molecule 3 (FAIM3) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60667
Sensitivity	0.058ng/mL
Entrez	9214
Expiration Date	Please enquire.