

Product Datasheet

Mouse FAM135B ELISA Kit (orb439446)

Description	Mouse Family With Sequence Similarity 135, Member B (FAM135B) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9DAI6
Sensitivity	0.057ng/mL
Entrez	70363
Expiration Date	Please enquire.