

Product Datasheet

Human FAM3B ELISA Kit (orb439448)

Description	Human Family With Sequence Similarity 3, Member B (FAM3B) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P58499
Sensitivity	0.235ng/mL
Entrez	54097
Expiration Date	Please enquire.