

Product Datasheet

Rat FAPa ELISA Kit (orb439453)

Description	Rat Fibroblast Activation Protein Alpha (FAPa) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8R492
Sensitivity	0.112ng/mL
Entrez	192203
Expiration Date	Please enquire.