

Product Datasheet

Mouse FAS ELISA Kit (orb439460)

Description	Mouse Factor Related Apoptosis (FAS) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P25446
Sensitivity	10.2pg/mL
Entrez	14102
Expiration Date	Please enquire.