

Product Datasheet

Monkey FAS ELISA Kit (orb439467)

Description	Monkey Factor Related Apoptosis (FAS) Sandwich ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9BDP2
Sensitivity	3.6pg/mL
Entrez	574332
Expiration Date	Please enquire.