

Product Datasheet

Rat FATP1 ELISA Kit (orb439468)

Description	Rat Fatty Acid Transport Protein 1 (FATP1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P97849
Sensitivity	0.071ng/mL
Entrez	94172
Expiration Date	Please enquire.