

Product Datasheet

Human FATP5 ELISA Kit (orb439469)

Description	Human Fatty Acid Transport Protein 5 (FATP5) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Y2P5
Sensitivity	0.055ng/mL
Entrez	10998
Expiration Date	Please enquire.