

Product Datasheet

Mouse FBLN2 ELISA Kit (orb439474)

Description	Mouse Fibulin 2 (FBLN2) Sandwich ELISA Kit
Reactivity	Mouse
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P37889
Sensitivity	21.45g/mL
Entrez	14115
Expiration Date	Please enquire.