

Product Datasheet

Human FcαR ELISA Kit (orb439484)

Description	Human Fc Fragment Of IgA Receptor (FcαR) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P24071
Sensitivity	0.35ng/mL
Entrez	2204
Expiration Date	Please enquire.