

Product Datasheet

Human FcεRII ELISA Kit (orb439485)

Description	Human Receptor II For The Fc Region Of Immunoglobulin E (FcεRII) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P06734
Sensitivity	0.107ng/mL
Entrez	2208
Expiration Date	Please enquire.