

Product Datasheet

Human FcgBP ELISA Kit (orb439488)

Description	Human Fc Fragment Of IgG Binding Protein (FcgBP) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y6R7
Sensitivity	0.116ng/mL
Entrez	8857
Expiration Date	Please enquire.