

Product Datasheet

Human FcγR3A ELISA Kit (orb439489)

Description	Human Fc Fragment Of IgG Low Affinity IIIa Receptor (FcγR3A) Sandwich ELISA Kit
Reactivity	Human
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P08637
Sensitivity	34ng/mL
Entrez	2214
Expiration Date	Please enquire.