

## Product Datasheet

### Human FcγR3B ELISA Kit (orb439490)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Fc Fragment Of IgG Low Affinity IIIb Receptor (FcγR3B) Sandwich ELISA Kit          |
| <b>Reactivity</b>        | Human                                                                                    |
| <b>Range</b>             | 1.25-80ng/mL                                                                             |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell lysates and other biological fluids.             |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates and other biological fluids.             |
| <b>Uniprot ID</b>        | <b>075015</b>                                                                            |
| <b>Sensitivity</b>       | 0.55ng/mL                                                                                |
| <b>Entrez</b>            | <b>2215</b>                                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                          |