

Product Datasheet

Rat FDPS ELISA Kit (orb439502)

Description	Rat Farnesyl Diphosphate Synthase (FDPS) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P05369
Sensitivity	0.057ng/mL
Entrez	83791
Expiration Date	Please enquire.