

Product Datasheet

Bovine FETUB ELISA Kit (orb439508)

Description	Bovine Fetuin B (FETUB) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q58D62
Sensitivity	11.7ng/mL
Entrez	504615
Expiration Date	Please enquire.