

Product Datasheet

Human FGa ELISA Kit (orb439514)

Description	Human Fibrinogen Alpha (FGa) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma and other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	P02671
Sensitivity	1.13ng/mL
Entrez	2243
Expiration Date	Please enquire.