

Product Datasheet

Human FGb ELISA Kit (orb439518)

Description	Human Fibrinogen Beta (FGb) Sandwich ELISA Kit
Reactivity	Human
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P02675
Sensitivity	29ng/mL
Entrez	2244
Expiration Date	Please enquire.