

Product Datasheet

Rat FGg ELISA Kit (orb439522)

Description	Rat Fibrinogen Gamma (FGg) Sandwich ELISA Kit
Reactivity	Rat
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P02680
Sensitivity	27ng/mL
Entrez	24367
Expiration Date	Please enquire.