

Product Datasheet

Human FGL1 ELISA Kit (orb439524)

Description	Human Fibrinogen Like Protein 1 (FGL1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q08830
Sensitivity	0.055ng/mL
Entrez	2267
Expiration Date	Please enquire.