

Product Datasheet

Mouse Flt3L ELISA Kit (orb439548)

Description	Mouse FMS Like Tyrosine Kinase 3 Ligand (Flt3L) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P49772
Sensitivity	6.9pg/mL
Entrez	14256
Expiration Date	Please enquire.