

Product Datasheet

Human FMO1 ELISA Kit (orb439549)

Description	Human Flavin Containing Monooxygenase 1 (FMO1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q01740
Sensitivity	0.115ng/mL
Entrez	2326
Expiration Date	Please enquire.