

Product Datasheet

Human FN3K ELISA Kit (orb439553)

Description	Human Fructosamine-3-Kinase (FN3K) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H479
Sensitivity	0.071ng/mL
Entrez	64122
Expiration Date	Please enquire.