

Product Datasheet

Mouse FOS ELISA Kit (orb439567)

Description	Mouse V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog (FOS) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P01101
Sensitivity	0.063ng/mL
Entrez	14281
Expiration Date	Please enquire.