

Product Datasheet

Rat FOS ELISA Kit (orb439568)

Description	Rat V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog (FOS) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P12841
Sensitivity	0.067ng/mL
Entrez	314322
Expiration Date	Please enquire.