

Product Datasheet

Rat FOXO1 ELISA Kit (orb439571)

Description	Rat Forkhead Box Protein O1 (FOXO1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	G3V7R4
Sensitivity	0.136ng/mL
Entrez	84482
Expiration Date	Please enquire.