

Product Datasheet

Human FOXO3 ELISA Kit (orb439572)

Description	Human Forkhead Box Protein O3 (FOXO3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O43524
Sensitivity	0.16ng/mL
Entrez	2309
Expiration Date	Please enquire.