

Product Datasheet

Mouse FOXO3 ELISA Kit (orb439573)

Description	Mouse Forkhead Box Protein O3 (FOXO3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9WVH4
Sensitivity	0.054ng/mL
Entrez	56484
Expiration Date	Please enquire.