

Product Datasheet

Human FPR1 ELISA Kit (orb439585)

Description	Human Formyl Peptide Receptor 1 (FPR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	cell lysates, cell culture supernates and other biological fluids.the Immunogen sequence for our orb439585 kit is C-terminal 267-350 AA.
Sample Types	cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P21462
Sensitivity	0.062ng/mL
Entrez	2357
Expiration Date	Please enquire.