

Product Datasheet

Porcine FSH ELISA Kit (orb439593)

Description	Porcine Follicle Stimulating Hormone (FSH) competitive inhibition ELISA Kit
Reactivity	Porcine
Range	6.25-400mIU/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01228
Sensitivity	2.13mIU/mL
Entrez	396895
Expiration Date	Please enquire.