

Product Datasheet

Rat FSH ELISA Kit (orb439594)

Description	Rat Follicle Stimulating Hormone (FSH) competitive inhibition ELISA Kit
Reactivity	Rat
Range	0.468-30mIU/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.21mIU
Expiration Date	Please enquire.