

Product Datasheet

Human FSIP1 ELISA Kit (orb439598)

Description	Human Fibrous Sheath Interacting Protein 1 (FSIP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8NA03
Sensitivity	0.112ng/mL
Entrez	161835
Expiration Date	Please enquire.