

Product Datasheet

Human FUCa2 ELISA Kit (orb439615)

Description	Human Fucosidase Alpha L2, Plasma (FUCa2) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9BTY2
Sensitivity	0.216ng/mL
Entrez	2519
Expiration Date	Please enquire.