

Product Datasheet

Human GABPa ELISA Kit (orb439637)

Description	Human GA Binding Protein Transcription Factor Alpha (GABPa) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q06546
Sensitivity	0.057ng/mL
Entrez	2551
Expiration Date	Please enquire.