

Product Datasheet

Human gABRa2 ELISA Kit (orb439638)

Description	Human Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P47869
Sensitivity	0.059ng/mL
Entrez	2555
Expiration Date	Please enquire.