

Product Datasheet

Mouse gABRa2 ELISA Kit (orb439639)

Description	Mouse Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P26048
Sensitivity	0.051ng/mL
Entrez	14395
Expiration Date	Please enquire.