

Product Datasheet

Rat GADD45a ELISA Kit (orb439645)

Description	Rat Growth Arrest And DNA Damage Inducible Protein Alpha (GADD45a) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q66HL6
Sensitivity	0.062ng/mL
Entrez	25112
Expiration Date	Please enquire.