

Product Datasheet

Mouse GADD45g ELISA Kit (orb439646)

Description	Mouse Growth Arrest And DNA Damage Inducible Protein Gamma (GADD45g) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z111
Sensitivity	0.225ng/mL
Entrez	23882
Expiration Date	Please enquire.