

Product Datasheet

Bovine GAL1 ELISA Kit (orb439653)

Description	Bovine Galectin 1 (GAL1) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P11116
Sensitivity	27.5pg/mL
Entrez	326598
Expiration Date	Please enquire.