

Product Datasheet

Rat GAL9 ELISA Kit (orb439682)

Description	Rat Galectin 9 (GAL9) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P97840
Sensitivity	35pg/mL
Entrez	25476
Expiration Date	Please enquire.