

Product Datasheet

Human GAP43 ELISA Kit (orb439687)

Description	Human Growth Associated Protein 43 (GAP43) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P17677
Sensitivity	0.062ng/mL
Entrez	2596
Expiration Date	Please enquire.