

Product Datasheet

Rat GAS6 ELISA Kit (orb439694)

Description	Rat Growth Arrest Specific Protein 6 (GAS6) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q63772
Sensitivity	0.071ng/mL
Entrez	58935
Expiration Date	Please enquire.