

Product Datasheet

Human GATA3 ELISA Kit (orb439695)

Description	Human GATA Binding Protein 3 (GATA3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P23771
Sensitivity	0.135ng/mL
Entrez	2625
Expiration Date	Please enquire.