

Product Datasheet

Human GbA ELISA Kit (orb439696)

Description	Human Glucosidase Beta, Acid (GbA) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P04062
Sensitivity	0.59ng/mL
Entrez	2629
Expiration Date	Please enquire.