

Product Datasheet

Human Glucagon ELISA Kit (orb439697)

Description	Human Glucagon (GC) competitive inhibition ELISA Kit
Reactivity	Human
Range	23.437-1500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01275
Sensitivity	7.3pg/mL
Entrez	2641
Expiration Date	Please enquire.