

Product Datasheet

Human GDA ELISA Kit (orb439711)

Description	Human Guanine Deaminase (GDA) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Y2T3
Sensitivity	0.102ng/mL
Entrez	9615
Expiration Date	Please enquire.