

Product Datasheet

Rat GDF1 ELISA Kit (orb439720)

Description	Rat Growth Differentiation Factor 1 (GDF1) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q1HI69
Sensitivity	12.5pg/mL
Entrez	306351
Expiration Date	Please enquire.